Mesh|EX[®] System & Sterilization Tray User Guide

EXSURCO MEDICAL
10804 GREEN ROAD
WAKEMAN, OH 44889

PN#: X101880, REV. A

Table of Contents

System Overview	1
Intended Use	2
Indications for Use.....	2
Description	2
Components	2
How Supplied.....	2
General Warnings and Precautions	2
Contraindications.....	3
Set Up Instructions.....	4
Inspection	4
Assembly	4
Operating Instructions	7
Skin and Tissue Graft Meshing	7
Postoperative Procedure	9
Cleaning and Sterilization Instructions	10
Manual Cleaning of the Mesh EX Reusable System	10
Automated Cleaning of the Mesh EX Reusable System.....	11
Lubrication	12
Sterilization	13
Troubleshooting	14
Warranty, Service, and Returns	14
Limited Warranty.....	15
Service	15
Returns.....	16
Appendix A: Technical Specifications	17
Appendix B: Parts Information.....	18
Appendix C: Glossary of Symbols.....	19
Help Information.....	20
Contact Information	20

System Overview

Intended Use

A reusable device that processes the skin from a whole piece of tissue with no fenestrations (holes) into a piece of tissue with apertures into the graft, allowing it to expand up to four times its size. This allows for greater surface area able to be covered with the graft. System may allow for changing meshing ratio or size of aperture in the graft, and max width of 4.25" graft. This may enable improved drainage, increased edge exposure, minimal contracture, and conformance to irregular body surfaces.

Indications for Use

A reusable manual device for meshing-skin (or tissue, substitute skin, allograft) in the operating room and clinical setting for use in wound, burn and other surgical procedures requiring meshed tissue for wound coverage including both inpatient and outpatient settings.

Description

The Exsurco Mesh|EX is a mechanically powered instrument capable of creating three unique mesh ratios on varying thickness of skin grafts from 0.005 to 0.020 inches (0.012-0.050 cm) and up to widths of 4.25" (10.80 cm). A separate cutting roller for each mesh ratio is included for use depending on tissue graft size needs.

The Mesh|EX Sterilization Tray will consist of an internal instrument basket with component placement bracketry, handles, and a lid. The brackets to hold the individual components are configured within the sterilization tray specifically for securing the Exsurco Mesher Components.

Components

1. Mesher Body
2. Feed Roller^a
3. Cutting Roller(s)
4. Graft Carrier^b
5. Sterilization Tray^c

a: also referred to as Propulsion Roller

b: Not Pictured

c: Filter sold separately

How Supplied

The system is supplied non-sterile and is reusable after cleaning and sterilization prior to use. All components must be inspected prior to cleaning and sterilization for any damages upon receipt.

Each piece of the system is shipped individually contained in protective boxes and must be opened, visually inspected, cleaned and then placed in the sterilization tray for sterilization.

The single-use, disposable skin graft carriers are provided sterile and sold separately. They should be visually inspected for packaging breaches or damage incurred during shipping.

WARNING: Re-sterilization or re-processing of the disposable skin graft carriers for reuse can result in serious physical harm to the patient.

NOTE: Only Exsurco approved parts & accessories are authorized for use with device. See Appendix B Parts Information.

General Warnings and Precautions

- The user of the Exsurco Mesher must be a trained professional and thoroughly familiar with the function, application, and instructions for use to avoid injury to the patient and operating room staff.
- Each physician must evaluate the appropriateness of the mesh ratio and technique for each patient based on their own medical training and expertise to prevent adverse patient outcomes. Each patient is different, and results can vary.
- Do not use cutting rollers or skin graft carriers from other manufacturers in the Exsurco Mesh|EX. Doing so would void the warranty of your product and could cause damage or harm to the user, device, and/or skin graft.
- Care should be taken when handling cutting rollers as blades are sharp and can cause injury.
- The operator of the Exsurco Mesh|EX and other clinical personnel should follow the universal protocol for aseptic technique as well as infection control as with all surgical procedures to prevent contamination to the device or patient.
- Do not use any other sterilization tray for the Exsurco Mesher system. If any other sterilization tray is used the sterilization of the Exsurco Mesh|EX cannot be guaranteed. The facility is responsible to complete an internal sterilization validation.
- Do not attempt to open the Exsurco Mesh|EX while the carrier is in the cutting mechanism, doing so may cause harm to the user, device, and skin graft.
- This product can expose you to the substance PTFE, which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov.

Contraindications

- The equipment should not be used for the preparation of tissue grafts > 1mm (0.04 inch) thickness.
- Meshed grafts may not be suitable for patients known to form keloids, or in areas subject to constant trauma. Use in these areas can result in severe scar contractures.

Set Up Instructions

Inspection

1. Ensure that sterilization is confirmed.
2. Verify all components as well as sterilization tray are free of damage and particulate. If any contaminants are present, re-clean the components.
3. Confirm cutting rollers are free from any bent/ chipped blade edges.
4. Verify all latches are secure, gears are engaging properly and are all in working order.
5. Inspect all movable parts of the system to ensure that they function through the intended range of motion.

Assembly

1. Remove the device from the sterilization tray and place on a sturdy flat surface. Recommended to be approximately waist height or 32" to 48" off the ground.
2. Place device in 'Open' position.
 - a. Pull the knobs on the latches out to release as depicted.


- b. Flip the top handle and cutting roller cage backwards to expose the roller groove.

3. Align rollers within roller groove.
 - a. Place bottom propulsion roller into groove.


- b. Set cutting roller on top, ensuring the gear on each roller are aligned correctly.

CAUTION: Cutting roller blades are sharp handle with care.

NOTE: Cutting rollers are color coded based on their expansion ratios. Utilize the below table to select the appropriate one. The rollers can also be identified by the laser etching on each gear.


Exsurco Mesh EX Cutting Roller	Expansion Ratio	Roller Color
Cutting Roller 'A'	1:1	Blue
Cutting Roller 'B'	2:1	Orange
Cutting Roller 'C'	3:1	Black

4. Close and lock the device.
 - a. Flip the top handle with cutting roller cage forward, encasing the rollers with the cage.


- b. Ensure the top is seated as far down as possible to align the latches and push the latches back into their locked position.

- c. Ensure the handle gear is aligned with the cutting roller gears and rotates as intended.

5. Open disposable graft carrier from sterile packaging.
 - a. Check the packaging integrity and expiration date. Ensure the carrier is free of particulate and foreign objects. Discard the graft carrier if it is damaged or defective.
 - b. Using the peel-pack end with peel corners, pull open the graft carrier pouch, allowing the sterile operator to reach in and grab the carrier.

Operating Instructions

Skin and Tissue Graft Meshing

- Each physician must evaluate the appropriateness of the technique for each patient based on their own medical training and expertise. Each patient is different, and results can vary.
 - It is recommended to begin meshing at low speed for better control of outcomes.
1. Place the tissue or recovered skin graft onto the graft carrier. For best results ensure that no parts of the graft are hanging off the carrier edges and that the graft is completely flat on the carrier.

NOTE: The tissue graft should be between 0.005" and 0.020" (0.012mm-0.050mm) thick.

2. Align the Mesh|EX Graft Carrier to the front of the Mesh|EX. Align one edge of the carrier flush to one side of the guide railings and move the graft carrier forward until the front contacts the cutting roller.

3. Place one hand on handle at top of device for stabilization during use.
4. Using the other hand, begin to rotate the crank handle in a clockwise direction.

NOTE: If the tissue graft is thicker, and therefore more difficult to process, it is advantageous to position the crank handle at the top of its rotation before allowing the graft carrier to contact the cutting roller. This allows the first motion to be a push rather than a pull which will make beginning the rotation easier.

NOTE: Do not attempt to open the Exsurco Mesh|EX while the graft carrier is in the cutting mechanism.

NOTE: Rotating the handle counterclockwise once the graft meshing has already started, may create misaligned cuts along the skin graft.

5. Continue rotation until graft and carrier are completely through and released by the cutting roller.


NOTE: Be aware of the potential for the tissue to 'stick' to the cutting roller. If this starts to happen, pause the rotation, and pull the graft back down onto the carrier before starting rotation back up again. Be careful when using metal instruments, as these can damage the cutting roller blades.

Postoperative Procedure

Perform this procedure after all the tissue grafts have been meshed for this operation. Dispose of all used carriers and sterilization tray filter in accordance with hospital policy for contaminated waste.

1. Open the Mesh|EX using the front latches and remove the cutting roller and propulsion roller.
2. Remove any apparent debris from the Mesher components and sterilization tray to prevent the drying of blood and body fluids.
3. Place the cutting rollers and the base into the sterilization container, in their appropriate location as depicted.

NOTE: It is important to be sure the cutting rollers are fully seated within their clamps to prevent movement and damage during transportation or sterilization.

4. Follow facility protocol to return the used Mesher and sterilization tray to Central Sterile Department (SPD) for cleaning and reprocessing.

CAUTION:

- Do not allow blood, tissue, or disinfectants to dry on any part of the Mesh|EX system.
- Do not use cleaning agents with chlorine or chloride (e.g., saline or bleach) as the active ingredient. They are corrosive to stainless steel and can result in physical damage to the system, including rust.
- Always use cleaning agents that are safe for surgical instruments and stainless steel.

Cleaning and Sterilization Instructions

The Mesh|EX system and sterilization tray are reusable and require cleaning and sterilization prior to reuse. Manual cleaning and automated cleaning have both been validated for use with the system. It is important to exercise care when handling the components to protect them from physical damage. Only remove the components as instructed in the following procedures. Any further disassembly of the system can cause damage and will void the warranty. The system contains no user serviceable parts and must be returned to Exsurco Medical for repair and maintenance.

As with any decontamination procedure, follow the accepted guidelines for hand washing, use of personal protective equipment (PPE) etc., as recommended by the Association of the Advancement of Medical Instrumentation (AAMI) Standards and Recommended practice, *Safe Handling and Biological Decontamination of Medical Devices in Health Care Facilities and in Non-Critical Settings*, current version of ANSI/AAMI ST35.

Always inspect instruments and devices for rust and corrosion per AORN and hospital procedures.

WARNING: Failure to comply to the below validated cleaning methods may result in corrosion of Mesh|EX components and could result in procedure delays or adverse outcomes. Please take care to follow instructions diligently and completely.

Manual Cleaning of the Mesh|EX Reusable System

1. Prepare an instrument-safe enzymatic cleaning solution which is compatible with aluminum and stainless steel according to the cleaning solution manufacturer specifications.
2. Immerse the Mesh|EX components in the cleaning solution and let soak for a minimum of fifteen (15) minutes. Do not immerse the unit overnight.
3. While still immersed in the cleaning solution, use a non-metal, soft-bristled brush moistened with cleaning solution to scrub each component until visibly clean.
4. Prepare a warm deionized water bath and immerse all components for a minimum of five (5) minutes. Do not immerse the unit overnight.
5. Rinse each component in warm deionized water for a minimum of one (1) minute.
6. Allow each component to fully dry. Utilizing pressurized air to blow liquid off cutting roller and speed up drying process is recommended. A lint-free wipe may be used on the propulsion roller or Mesh|EX body if desired.
7. In a well-lit area, visually inspect each component to ensure all surfaces are clean and dry.
8. Lubricate the device following the below instructions.

Automated Cleaning of the Mesh|EX Reusable System

If automated cleaning is available at the facility complete the cleaning per the automated cycle parameters validated below:

1. Prepare an instrument-safe enzymatic cleaning solution which is compatible with aluminum and stainless steel according to the cleaning solution manufacturer specifications.
2. Immerse the Mesh|EX components in the cleaning solution and let soak for a minimum of fifteen (15) minutes. Do not immerse the unit overnight.
3. While still immersed in the cleaning solution, use a non-metal, soft-bristled brush moistened with cleaning solution to scrub each component until visibly clean.
4. Rinse each component in deionized water for a minimum of one minute and 30 seconds (1:30).
5. Place the disassembled components into a wire basket, excluding the sterilization tray components which should be placed directly into the chamber.

NOTE: The Sterilization Tray base should be placed with the cavity facing downward as to not accumulate water.

Table 1 (Automated Cleaning Instructions for Reusable System)

Step	Time (Minutes)	Recommended Temperature
Pre-wash	3:00	60°C
Cleaning ^a	3:00	85°C
Neutralizing ^b	1:00	20°C
Rinse I	1:30	20°C
Rinse II	1:30	20°C
Final Rinse ^c	2:00	80°C
Drying	25:00	110°C

a: Washer cleaning solution should be an enzymatic solution compatible with aluminum.

b: Neutralizing solution should be appropriate for the utilized cleaning solution, based upon the manufacturers' recommendation.

c: Final Rinse should be completed with critical water.

Upon completion of the automated cycle, inspect as outlined below:

1. In a well-lit area, visually inspect each component to ensure all surfaces are clean and dry.
2. If any contaminants are present, re-clean following the appropriate cleaning instructions. Check for any corrosion on the laser etching of the components. If present, contact Exsurco Medical Customer Service for a Service Repair Order at 1-800-243-6049 ext. 4, or email support@exsurco.com.
3. Lubricate the device following the below instructions.


Lubrication

It is important to add lubrication to the device to ensure latches and hinges continue to operate freely and easily over time.

1. Use a water-soluble surgical grade instrument lubricant (such as instrument milk).

WARNING: Silicone-based or oil-based lubricants should never be used during sterilization as they can cause interference with steam and prevent effective sterilization.

2. After cleaning and before sterilization, spray or submerge the Mesh|EX components in lubrication. Ensure lubrication is applied to the cutting roller blades, hinges, bearings, latches, pinion and other moving parts as pictured below.

3. Sterilize per the below specifications.

Sterilization

The Exsurco Mesh|EX system has been validated for steam sterilization using the Exsurco Mesh|EX sterilization tray. Should any other method or sterilization tray be used, it is the sole responsibility of the facility to validate any sterilization parameters that are not provided directly by the manufacturer.

Complete cleaning process as directed above before beginning sterilization.

WARNING: Failure to follow to procedures outlined in this section may result in non-sterile product. The facility is responsible for ensuring all their personnel are familiar with the cleaning and sterilization procedures.

1. Secure the cleaned Mesh|EX system components into their proper locations in the cleaned sterilization tray.
2. Utilize facility approved and required sterilization indicators.
3. Replace the filter contained within the lid of the sterilization tray.

NOTE: Use facility approved filters, Exsurco Medial does not provide filter replacements. Do not use damaged, torn, or wet filters. The filter within the sterilization tray is disposable and should be replaced after each use. The following are suitable for usage with the sterilization tray; Paper, Textile, PTFE (Teflon), or PP (polypropylene) filters.

4. Place the lid on top of the base, and then secure the side latches.

The Mesh|EX system is ready for sterilization. Only use an autoclave validated in compliance with current version of AAMI ST8, Hospital Steam Sterilizers.

Utilize the validated steam sterilization method listed below to provide a sterility assurance level of 10^{-6} .

Cycle Type	Pre-Vacuum
Set-Point Temperature	270°F
Exposure Time	4 minutes
Drying Time	30 minutes

NOTE: Pre-Vacuum Steam, stacking should not exceed 18 inches in height for effective air removal and adequate steam penetration.

NOTE: Temperature, Exposure and Drying Times listed above are minimum values which are required to achieve a 10^{-6} sterility assurance level (SAL).

Upon completing sterilization, store the Mesh|EX system in the sterilization tray in a dry and dust-free place. The shelf life is dependent on the sterile barrier employed, manner in which product is stored, environment, and handling conditions. A maximum shelf life for sterilized medical devices before use should be defined by each healthcare facility.

NOTE: References & Internal Validation details are on file at Exsurco Medical and available upon request.


Troubleshooting

Use the following assessment steps to confirm proper operation of the Exsurco Mesh|EX and Sterilization Tray. If you discover an issue with any component, do not use the system or sterilization tray and return it to Exsurco Medical for service. For additional information refer to the Warranty, Service and Returns section.

If you are experiencing a problem the first step should be to inspect the device, look for any of the following:

- Signs of physical damage or wear such as loose parts, missing or loose screws, denting or signs of being dropped, areas of corrosion.
- Hinges and gears continue to have their full range of mobility.
- Debris blocking any areas, preventing component interaction.

If no issues are observed during this assessment yet problem persists, see the troubleshooting table below:

Problem	Possible Solution
Operating the Mesh EX is strenuous	The Mesh EX Graft Carrier may not be aligned properly within the guard rails. Check alignment before proceeding.
The cutting roller does not rotate, while turning the handle	<p>The gears may not be engaged correctly. Ensure that the top of the mesher is completely seated and the latches are in a fully closed position. Verify the cutting roller gears and the propulsion roller gears are mated correctly.</p> 
The graft carrier is not moved through the Mesh EX, while turning the handle	Verify that the skin graft carrier is flush with the cutting roller and is engaging with the blade edges.
The skin graft is only cut partially or not at all	Cutting roller blades may be damaged or dull. Inspect the cutting roller for any defects. If blades are damaged, replace the cutting roller.
Graft migrates on carrier during operation	Check the thickness of the skin graft. Subcutaneous tissue included on the skin graft can cause adverse results. Skin graft may be too saturated with saline solution, attempt to remove excess liquid.
Deformation or rupture of the internal lid gasket occurred	<p>Replacement of the gasket needs to be completed per manufacturers usage instructions:</p> <ol style="list-style-type: none"> 1. Clean the lid groove with alcohol to remove any gasket residues. 2. Apply silicone adhesive (Henkel-Sista® recommended) to the gasket groove. 3. Cut a 136cm gasket and place it into the groove (flat surface facing up) and wait 12 hours to completely dry.
The handle gear is aligned with the cutting roller gears and does not rotate as intended	Inspect the gear rollers for damage that may have occurred during reprocessing. Ensure nothing is jammed within the pinion or between any gears preventing proper rotation.

Warranty, Service, and Returns

This section describes the limited warranty, service requirements, and returns policy.

Limited Warranty

Exsurco Medical (“Exsurco”) warrants that each new Mesh|EX system, sterilization tray, and its components have been tested and inspected and have left the facility in proper working condition. Exsurco warrants the product is free from manufacturing defects in material and workmanship under normal use and service for a period of one year from the date of delivery by Exsurco to the customer. The sterile, disposable graft carrier is not warranted beyond the expiration date stated on the product labeling.

The Exsurco mesher and sterilization tray are guaranteed to be free from defects in material and workmanship when maintained by the customer, properly cleaned, and used under normal circumstances for its intended purpose.

During the warranty period, Exsurco shall repair, or at its sole opinion, replace the defective product or part without cost incurred to the customer. Defective parts replaced under this warranty shall become the property of Exsurco. This warranty does not cover damage caused by misuse, abuse (dropped or damaged), accident, neglect, or any use not prescribed in this user guide. If the system becomes defective because of misuse or abnormal conditions of operations, these repairs shall be charged to the customer. All products returned for warranty or non-warranty repair must be through the Exsurco return procedure. For more information, refer to the Returns information below. In no event shall Exsurco be liable for any incidental, indirect, consequential, or punitive damages in connection with the acquisition or use of Exsurco product. Further, this warranty shall not apply to, and Exsurco shall not be responsible for, any loss or damage arising in connection with the purchase or use of Exsurco product that has been repaired by anyone other than Exsurco, or altered in any way so as to, in the judgement of Exsurco Medical, affect the usability or reliability of the product. If this happens, the product repairs shall be charged to the customer. This limited warranty is exclusive and in lieu of all other warranties, expressed or implied, and of all other obligations or liabilities on Exsurco’s part, and Exsurco neither assumes nor authorizes any representative or other person to assume for it any other liability in connection with this product.

Exsurco disclaims all other warranties, expressed or implied, including any implied warranty of merchantability or of fitness for a particular purpose or application or warranty of quality as well as any expressed or implied warranty to patients.

Service

The Exsurco Mesh|EX system and sterilization tray contains no serviceable parts. The system must be returned to Exsurco Medical for servicing. Exsurco Medical trained Technical Service Technicians perform all repairs in-house.

NOTE: Exsurco Medical cannot be held liable for any system malfunctions resulting from repairs or service performed by anyone other than Exsurco Medical.

Do not attempt to open, remove any screws, or disassemble any components. Tampering with hardware will result in altered calibration and will subsequently affect cutting performance. If this happens the warranty is immediately voided, and the customer must return the system to Exsurco Medical for repair that shall be charged to the customer.

Our Technical Service Department can be contacted at 1-800-243-6049, ext. 5 or email support@exsurco.com to answer any general questions that relate to the operation, warranty, and service of the system.

The Exsurco Mesh|EX system must be returned annually for inspection and preventative maintenance. Annual factory maintenance is strongly recommended in order to ensure consistent operation and accuracy.

Returns

When it is necessary to return the Exsurco Mesh|EX or sterilization tray for inspection, preventative maintenance, warranty service, or non-warranty service, please contact Customer Service at 1-800-243-6049, ext. 4, or email support@exsurco.com to receive a Return Material Authorization (RMA) or Service Repair Order (SRO).

Product will not be accepted without an RMA/SRO. Exsurco Medical can provide the customer with the approved return shipping container. The customer is responsible for the return freight.

When the customer is returning the system, we require proof of sterilization documents upon receipt of the RMA or SRO. If proof of sterilization is not included with the system, the system shall be returned for sterilization. Exsurco Medical does not accept product that has been used without being sterilized first or has not been shipped with the sterilization documents. If it is a new product or has not been used, then customer service shall note that on the RMA/SRO documentation.

Upon receiving the returned product at Exsurco, the product is reviewed by our Technical Service Department. Customer Service shall provide the customer with a service estimate for any repairs not covered under warranty or extended warranty. Product will be returned to the customer in approved protective shipping container.

Exsurco Medical stands behind the quality of our products. Should you have a complaint or not be completely satisfied with your product, please call our Customer Service team. They are available Monday through Friday, 8am to 5pm EST and is here to assist you in any way. Email support@exsurco.com or call 1-800-243-6049, ext. 4. To help investigate the event further, we require the following information:

- Name, address, and account number.
- Phone number and email address.
- Part/product that the complaint is being filed on.
- Serial number (equipment) or lot number (disposable).
- A description of the event. Be as descriptive as possible, and if available, provide pictures.

Once a complaint is received, we may contact you for more information. When the product is returned, a preliminary evaluation is generally conducted, and the product is repaired (if required) and restored to conformance prior to being shipped back. A closure letter shall be provided advising you of the outcome of the investigation regarding the reported event.

Appendix A: Technical Specifications

Exsurco Mesh|EX System

Specification	Value
Storage and Transport	As per internal hospital procedure.
Weight and Measurement	Assembled Mesh EX Weight: 8.34 lbs. (3.78 kg) Complete Mesh EX system in sterilization tray: 19.82 lbs. (8.99 kg) Assembled Mesh EX Dimensions: 11.5 x 5.5 x 6.5 in (29.21 x 13.97 x 16.51 cm)
Sterile, Single-Use Graft Carriers Environmental Conditions	It is recommended these materials be stored in the provided packaging in a cool (under 32°C (90°F)), clean and dry area. Humidity Condition Storage Range: 30-70%
Mesh EX System Operation Environmental Conditions	Indoors under normal Operating Room conditions with ambient lighting or better.
Skin Graft Parameters	Maximum Skin Graft Thickness: 0.020 inches (0.050 cm) Maximum Skin Graft Width: 4.25" (10.80 cm) Maximum Skin Graft Length: 13.8" (35.0 cm)
Mesh EX Meshing Ratios	A Cutting Roller: 1:1 Expansion Ratio B Cutting Roller: 2:1 Expansion Ratio C Cutting Roller: 3:1 Expansion Ratio
Disposal	If the Cutting Rollers, Propulsion Roller, or any part of the Mesh EX system is being discarded, it cannot be discarded into trash. Contact Exsurco Medical for assistance in this process.

Exsurco Mesh|EX Sterilization Tray

Specification	Value
Sterilization	Steam sterilization only in validated Mesh EX Sterilization Tray.
Cycle	Pre-vacuum steam sterilization
Storage and Transport	As per internal hospital procedure.
Weight	Lid: 2.28 lbs. (1.03 kg) Base: 4.26 lbs. (1.93 kg) Internal Instrument Basket: 1.92 lbs. (0.87 kg)
Dimensions	19 x 11 x 5.5 in (48.26 x 27.94 x 13.97 cm)
Disposal	If the Lid, Base, or Basket is being discarded, it cannot be discarded into trash. Contact Exsurco Medical for assistance in this process.

Appendix B: Parts Information

Part #: X101892

Description: Exsurco Mesher Kit

Quantity	Part Number	Description
1 each	X101933	Mesh EX Base
1 each	X101836	Mesh EX Body (included in X101933)
1 each	X101876	Mesh EX Feed Roller (included in X101933)
1 each	X101855	Mesh EX Cutting Roller, A
1 each	X101863	Mesh EX Cutting Roller, B
1 each	X101866	Mesh EX Cutting Roller, C
1 each	X101882	Mesh EX Sterilization Tray

Sterile, Disposable Graft Carriers sold separately.

NOTE: Use only the Skin Graft Carriers below for this product. Failure to do so voids the warranty and can cause injury.

Part Number	Description
X101886	Mesh EX Graft Carrier, 14" (20/box)

Ordering Information

To place an order for the Exsurco Mesh|EX system, please contact Customer Service at 1-800-243-6049, ext. 4, or email support@exsurco.com. For a price quote and all product-related questions, please contact your local Exsurco Medical account manager or our sales department at 1-800-243-6049, ext. 1, or email sales@exsurco.com.

The Exsurco Mesh|EX components are only sold as a replacement product to current customers. To order replacement components separately, please contact Customer Service at 1-800-243-6049, ext. 4, or email support@exsurco.com.

Appendix C: Glossary of Symbols

Symbol	Definition
	Manufacturer.
	Use by date. This symbol is intended to indicate that the system should not be used after the end of the date shown.
	Do not reuse. Single use, use only once.
	Symbol for method of sterilization using radiation.
	Consult instruction for use.
	Weight. Used to indicate mass.
	Do not use if package is damaged.
	Keep product dry.
	Caution: USA federal law restricts this system to sale by or on the order of a physician.
	Caution.
	Not made with natural rubber Latex.
	Warning. Sharp Element.
	Non-Sterile.

Help Information

Contact Information

Please refer to the information below for contacting Exsurco Medical directly for any of your questions, comments, or concerns.

Phone: 1-800-243-6049

Website: www.exsurco.com

Press the number below for the following departments or email at:

- 1- Sales & Clinical Support: sales@exsurco.com
- 2- Marketing & Product Inquiries: information@exsurco.com
- 4- Customer Service or to Report a Quality Concern: support@exsurco.com
- 5- System or Technical Service: support@exsurco.com
- 6- Quality Department: support@exsurco.com
- 7- Accounting Department: support@exsurco.com

Exsurco Medical
10804 Green Road
Wakeman, OH 44889

Exsurco Medical
P.O. Box 376
Vermillion, OH 44089

2023 Exsurco Medical, Inc. All Rights Reserved
Exsurco® and Mesh|EX® are registered trademarks of Exsurco Medical, Inc.
PN: X101880 Rev. A
GTIN: (01) 0 081000678013 0

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2023 Exsurco Medical, Inc. Tous droits réservés.
Exsurco® et Mesh|EX® sont des marques de commerce déposées d'Exsurco Medical, Inc.