


PREMIUM  
PRECISE  
EFFICIENT



**DermaSharp™ Replacement Blades with Edge-Retention Technology™ for longer-lasting sharpness**


**30%**  
Sharper than standard dermatome blades

Designed specifically with tissue banks in mind, Exsurco® Medical DermaSharp dermatome blades help facilities maximize efficiencies during allograft recoveries.

DermaSharp blades are ground and honed to a finer, sharper edge than standard dermatome blades. Then a unique coating process – Edge-Retention Technology – is applied, which preserves the blade edge during recoveries, resulting in fewer blade changes.


**Initial Sharpness**

Even straight out of the package, DermaSharp blades are 30% sharper than standard dermatome blades. See chart below.



**Edge Retention**

After 20 test cycles using abrasion resistant materials, DermaSharp blades remain 3 times sharper than standard dermatome blades. See chart below.



\*Sharpness numeric values are used for relative comparison.

**Cadaveric tissue bank use only**

Test documentation on file.



# EXSURCO® DERMASHARP™

## REPLACEMENT BLADES

- 30% sharper than standard dermatome blades
- Remain 3 times sharper than standard dermatome blades
- Feature Edge-Retention Technology™
- Precision-engineered for accurate control when recovering uniform grafts
- Premium quality, hardened, surgical stainless steel blades for efficient and consistent recovery
- Manufactured to exact specifications for use in popular competitor dermatome models
- Sterile-packed for tissue recovery settings
- Securely packaged for safe operator handling
- Sold in boxes of 10

Order today at [shop.exsurco.com](http://shop.exsurco.com)



Exsurco Part #X100448

### Exsurco Part #X100448

Length: 4 inches

Suitable for use only with PADGETT® Tissue Bank Dermatome:

- 3539-800 (Electric Model SB)



Exsurco Part #X100449

### Exsurco Part #X100449

Length: 4 inches

Suitable for use only with ZIMMER® Tissue Bank Dermatome:

- 00-8831-001-00 (Electric Model)

**Cadaveric tissue bank use only**



## The Edge-Retention Technology™ Difference

Extensive Exsurco Medical testing has shown DermaSharp Blades to be 30% sharper right out of the package and remain 3 times sharper than standard dermatome blades after use. Our Edge-Retention Technology - a unique combination of materials and processing - preserves the blade edge during recoveries, resulting in fewer blade changes.

Test documentation on file.

## Contact Us:

### Exsurco Medical, Inc.

6801 State Route 60  
Birmingham, OH 44889  
Phone: 800-243-6049  
Fax: 440-328-4535

[www.exsurco.com/dermasharp](http://www.exsurco.com/dermasharp)  
Sales: [sales@exsurco.com](mailto:sales@exsurco.com)  
Information: [info@exsurco.com](mailto:info@exsurco.com)  
Service and Support: [support@exsurco.com](mailto:support@exsurco.com)

EM0005 Rev. 1

Edge-Retention Technology is a trademark of Exsurco Medical, Inc.  
Exsurco is a registered trademark of Exsurco Medical, Inc.

Padgett is a registered trademark of Integra LifeSciences Corporation.  
Zimmer is a registered trademark of Zimmer, Inc.

©2013 Exsurco Medical. All rights reserved.

**exsurco**<sup>®</sup>  
medical