

**Part Numbers: 4Med Mesher and Accessories****For US Tissue Banks and Hospitals Only**

Vendor P/N	Description	Details	Carriers Used
<b>Reusable Mesher Product</b>			
9910100	Skin Graft Mesher with Ratchet, Variable Ratio	Variable Ratio (1:1 to 1:5)	Standard, Long
9910100-16	Skin Graft Mesher with Ratchet, Close Fixed Ratio 1:1	Close Fixed Ratio 1:1	Standard, Long
9910100-17	Skin Graft Mesher with Ratchet, 17cm, Close Fixed Ratio 1:1 (Wide)	Close Fixed Ratio 1:1 (Wide)	Wide, Standard, Long
Special Order	Skin Graft Mesher with Ratchet, Fixed 1:2	Call Exsurco for Details	Standard, Long
Special Order	Skin Graft Mesher with Ratchet, Fixed 1:3	Call Exsurco for Details	Standard, Long
Special Order	Skin Graft Mesher with Ratchet, Fixed 1:4	Call Exsurco for Details	Standard, Long
Special Order	Skin Graft Mesher with Ratchet, Fixed 1:5	Call Exsurco for Details	Standard, Long
<b>Disposable- Sterile Carriers</b>			
9910102	Skin Graft Carrier, Standard	Carrier Size (W x L): 14cm x 30cm (5.5" x 12.0")	N/A
9910102-60	Skin Graft Carrier, Long	Carrier Size (W x L): 14cm x 60cm (5.5" x 24.0")	N/A
9910102-35	Skin Graft Carrier, Wide	Carrier Size (W x L): 16.5cm x 35.5cm (6.6" x 14.2")	N/A
<b>Accessories</b>			
9910104	Sterilizing Container with Top Filter	Sterilization Container	N/A
9911010	Ratchet for Mesher	Ratchet	N/A
Y800.00	Filter for Sterilization Container	filters that fit the Sterilizing Container	N/A