



## Certificate of Conformance

Certificate Number: EXMJ2387X87

Dated Issued: 04/02/2019

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

<b>Product Item Number:</b>	X100457
<b>Item Description:</b>	ASSY, BOX OF 5 AMALGATOME MD PKG
<b>Item Lot Number:</b>	X-0000000000087
<b>Expiration Date:</b>	09/2019

The blades were manufactured using a validated sterilization process:

<b>STERIS Gamma Process Run ID:</b>	103819A		
<b>Processing Run Date:</b>	31-MAR-2019		
<b>Minimum Specified Dose (kGy):</b>	25.0	<b>Minimum Delivered Dose (kGy):</b>	28.7
<b>Maximum Specified Dose (kGy):</b>	45.0	<b>Maximum Delivered Dose (kGy):</b>	36.5

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**

Rushell Ringle  
Quality Department

04-03-2019  
Date

Swasita Saigal ON  
Quality Manager *behalf of Kayli Camp*

04-03-2019  
Date

<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.