



Certificate of Conformance

Certificate Number: EXMJ2549X92

Dated Issued: 24/JUL/2019

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	ASSY, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-0000000000092
Expiration Date:	06-2021

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	105731A		
Processing Run Date:	23-JUL-2019		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	27.7
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	36.2

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**



 Quality Department

24 Jul 2019

 Date



 Quality Manager

24 Jul 2019

 Date

¹ Manufacturing data filed within the Quality Management System and available onsite.