



Certificate of Conformance

Certificate Number: EXMXJ2501X34

Dated Issued: 08-JUL-2019


The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

| | |
|-----------------------------|--------------------------------|
| Product Item Number: | X100593 |
| Item Description: | Pack, Amalgatome SK head, 5 pc |
| Item Lot Number: | X-00000000000034 |
| Expiration Date: | MAY-2021 |

The blades were manufactured using a validated sterilization process:

| | | | |
|--------------------------------------|-------------|--------------------------------------|------|
| STERIS Gamma Process Run ID: | 105371A | | |
| Processing Run Date: | 28-JUN-2019 | | |
| Minimum Specified Dose (kGy): | 25.0 | Minimum Delivered Dose (kGy): | 28.3 |
| Maximum Specified Dose (kGy): | 45.0 | Maximum Delivered Dose (kGy): | 37.8 |

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.



 Quality Department

08-Jul-2019

 Date



 Quality Manager

08-July-2019

 Date

¹ Manufacturing data filed within the Quality Management System and available onsite.