



## Certificate of Conformance

Certificate Number: EXMXJ2660X94

Dated Issued: 10/01/2019

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

<b>Product Item Number:</b>	X100457
<b>Item Description:</b>	ASSY, BOX OF 5 AMALGATOME MD PKG
<b>Item Lot Number:</b>	X-0000000000094
<b>Expiration Date:</b>	2021/08

The blades were manufactured using a validated sterilization process:

<b>STERIS Gamma Process Run ID:</b>	106837A		
<b>Processing Run Date:</b>	29-SEPT-2019		
<b>Minimum Specified Dose (kGy):</b>	25.0	<b>Minimum Delivered Dose (kGy):</b>	27.5
<b>Maximum Specified Dose (kGy):</b>	45.0	<b>Maximum Delivered Dose (kGy):</b>	38.5

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.

  
Quality Department

10-01-2019  
Date

  
Quality Manager

10-03-2019  
Date

<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.