



Certificate of Conformance

Certificate Number: EXMXJ2662X39

Dated Issued: 10/01/2019

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100593
Item Description:	Pack, Amalgatome SK head, 5 pc
Item Lot Number:	X-00000000000039
Expiration Date:	2021/08

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	106837A		
Processing Run Date:	29-SEPT-2019		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	27.5
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	38.5

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.


Quality Department

10-01-2019
Date


Quality Manager

10-04-2019
Date

¹ Manufacturing data filed within the Quality Management System and available onsite.