



Certificate of Conformance

Certificate Number: EXMXJ3083X105

Dated Issued: 07/07/2020

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	ASSY, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-0000000000105
Expiration Date:	2022/05

The blades were manufactured using a validated sterilization process:


STERIS Gamma Process Run ID:	110726A		
Processing Run Date:	27-Jun-2020		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	27.7
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	38.1

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company.
Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**


 Quality Department

 Quality Manager


 Date


 Date

¹ Manufacturing data filed within the Quality Management System and available onsite.