



Certificate of Conformance

Certificate Number: EXMJ3364X112

Dated Issued: 12/30/2020

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	ASSY, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-0000000000112
Expiration Date:	2022/11

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	112659A		
Processing Run Date:	28-Dec-2020		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	28.2
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	40.2

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.


Quality Department


Quality Manager
for Rayli Camp

12/30/2020
Date

12-30-2020
Date

¹ Manufacturing data filed within the Quality Management System and available onsite.