



## Certificate of Conformance

Certificate Number: EXMJ3401X75

Dated Issued: 03/17/2021

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

<b>Product Item Number:</b>	X100593
<b>Item Description:</b>	Pack, Amalgatome SK head, 5 pc
<b>Item Lot Number:</b>	X-00000000000075
<b>Expiration Date:</b>	2023/02


The blades were manufactured using a validated sterilization process:

<b>STERIS Gamma Process Run ID:</b>	1073-2594A		
<b>Processing Run Date:</b>	08-Mar-2021		
<b>Minimum Specified Dose (kGy):</b>	25.0	<b>Minimum Delivered Dose (kGy):</b>	27.5
<b>Maximum Specified Dose (kGy):</b>	45.0	<b>Maximum Delivered Dose (kGy):</b>	39.3


Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**

  
 \_\_\_\_\_  
 Quality Department

3/17/21  
 \_\_\_\_\_  
 Date

  
 \_\_\_\_\_  
 Quality Manager

03-18-2021  
 \_\_\_\_\_  
 Date

Sara Ann Mackinlay  


<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.