



Certificate of Conformance

Certificate Number: EXMXJ3460X81

Dated Issued: 05/11/2021


The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100593
Item Description:	Pack, Amalgatome SK Head, 5 PC
Item Lot Number:	X-00000000000081
Expiration Date:	2023-04

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	1073-3685A		
Processing Run Date:	10-May-21		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	27.2
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	39.8


Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**



 Quality Department

05/13/2021

 Date



 Quality Manager

05-11-2021

 Date

¹ Manufacturing data filed within the Quality Management System and available onsite.