



Certificate of Conformance

Certificate Number: EXMXJ3689X125

Dated Issued: 11/12/2021

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	BLADE, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-00000000000125
Expiration Date:	2023/09

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	1073-6325A		
Processing Run Date:	01-Nov-2021		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	28.2
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	38.7

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company.

Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.



 Quality Department

11/12/2021

 Date



 Quality Manager

11/16/21

 Date

¹ Manufacturing data filed within the Quality Management System and available onsite.