



## Certificate of Conformance

Certificate Number: EXMXJ4199X140

Dated Issued: 01/19/2023

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

|                             |                                    |
|-----------------------------|------------------------------------|
| <b>Product Item Number:</b> | X100457                            |
| <b>Item Description:</b>    | BLADE, BOX OF 5, AMALGATOME MD PKG |
| <b>Item Lot Number:</b>     | X-0000000000140                    |
| <b>Expiration Date:</b>     | 2024/11                            |

The blades were manufactured using a validated sterilization process:

|                                      |             |                                      |      |
|--------------------------------------|-------------|--------------------------------------|------|
| <b>STERIS Gamma Process Run ID:</b>  | 1073-12322A |                                      |      |
| <b>Processing Run Date:</b>          | 02-Jan-2023 |                                      |      |
| <b>Minimum Specified Dose (kGy):</b> | 25.0        | <b>Minimum Delivered Dose (kGy):</b> | 28.6 |
| <b>Maximum Specified Dose (kGy):</b> | 45.0        | <b>Maximum Delivered Dose (kGy):</b> | 38.4 |


Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**



Quality Department



Quality Manager



Date



Date

<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.