



## Certificate of Conformance

Certificate Number: EXMXJ4695X152

Dated Issued: 2/8/2024

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

|                             |                                   |
|-----------------------------|-----------------------------------|
| <b>Product Item Number:</b> | X100457                           |
| <b>Item Description:</b>    | BLADE, BOX OF 5 AMALGATOME MD PKG |
| <b>Item Lot Number:</b>     | X-0000000000152                   |
| <b>Expiration Date:</b>     | 2026-01                           |

The blades were manufactured using a validated sterilization process:

|                                      |             |                                      |      |
|--------------------------------------|-------------|--------------------------------------|------|
| <b>STERIS Gamma Process Run ID:</b>  | 1073-17451A |                                      |      |
| <b>Processing Run Date:</b>          | 07-Feb-2024 |                                      |      |
| <b>Minimum Specified Dose (kGy):</b> | 25.0        | <b>Minimum Delivered Dose (kGy):</b> | 27.6 |
| <b>Maximum Specified Dose (kGy):</b> | 45.0        | <b>Maximum Delivered Dose (kGy):</b> | 38.4 |

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**

  
\_\_\_\_\_  
Quality Department

  
\_\_\_\_\_  
Quality Manager

2/9/24  
Date

02-15-2024  
Date

<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.