



Certificate of Conformance

Certificate Number: EXMJ5013X154

Dated Issued: 11/25/2024

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	ASSY, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-0000000000154
Expiration Date:	10/2026

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	1073-21221A		
Processing Run Date:	11/22/2024		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	27.5
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	39.1

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**

Heleen Con

Quality Department

Abdulhadi Alshaykh

Quality Manager

11/25/2024

Date

11/25/2024

Date

¹ Manufacturing data filed within the Quality Management System and available onsite.