



Certificate of Conformance

Certificate Number: EXMXJ5267X155

Dated Issued: 03/19/2025

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100457
Item Description:	ASSY, BOX OF 5 AMALGATOME MD PKG
Item Lot Number:	X-0000000000155
Expiration Date:	02-2027

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	1073-22774A		
Processing Run Date:	03/18/2025		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	28.8
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	39.0

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.

Rebecca Cove
Quality Department

[Signature]
Quality Manager

03/19/25
Date

3/19/25
Date

¹ Manufacturing data filed within the Quality Management System and available onsite.