



## Certificate of Conformance

Certificate Number: EXMJ5576X134

Dated Issued: 10/29/2025

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.<sup>1</sup>

<b>Product Item Number:</b>	X101835
<b>Item Description:</b>	ASK, CUTTING HEAD, 10 PC
<b>Item Lot Number:</b>	X-0000000000134
<b>Expiration Date:</b>	2027-09

The blades were manufactured using a validated sterilization process:

<b>STERIS Gamma Process Run ID:</b>	1073-25815A		
<b>Processing Run Date:</b>	10/27/2025		
<b>Minimum Specified Dose (kGy):</b>	25.0	<b>Minimum Delivered Dose (kGy):</b>	27.9
<b>Maximum Specified Dose (kGy):</b>	45.0	<b>Maximum Delivered Dose (kGy):</b>	36.9

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company.  
Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.

Rebecca Cover  
Quality Department

James E. Tanager  
Quality Manager

10/30/25  
Date

10/30/25  
Date

<sup>1</sup> Manufacturing data filed within the Quality Management System and available onsite.