



Certificate of Conformance

Certificate Number: EXMXJ5633X62

Dated Issued: 12/11/2025

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

| | |
|-----------------------------|------------------------------------|
| Product Item Number: | X101788 |
| Item Description: | PACK, AMALGATOME SD 4" BLADE, 5 PC |
| Item Lot Number: | X-0000000000062 |
| Expiration Date: | 11/17/2025 |

The blades were manufactured using a validated sterilization process:

| | | | |
|--------------------------------------|-------------|--------------------------------------|------|
| STERIS Gamma Process Run ID: | 1073-26317A | | |
| Processing Run Date: | 12/06/2025 | | |
| Minimum Specified Dose (kGy): | 25.0 | Minimum Delivered Dose (kGy): | 28.9 |
| Maximum Specified Dose (kGy): | 45.0 | Maximum Delivered Dose (kGy): | 38.4 |

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE.**

Rebecca Cover
Quality Department

12/11/25
Date

James E. Thompson
Quality Manager

12/11/25
Date

¹ Manufacturing data filed within the Quality Management System and available onsite.